



BEESWIFT
FOCUSED ON SAFETY

Product Code

CM0565

Available in Sizes

70ML

Available in Colours



White
CM0565

DATASHEET

Click Medical Cuteeze Haemostatic Spray 70ML



First aid training and regulations now include the use of haemostatic dressings and tourniquets for workers who are in high risk environments where serious trauma injuries and bleeding are a possibility.

Cut-eeze is our new range of haemostatic dressings

Description : Our haemostatic(or wound clotting) dressings comprise of a white knitted fabric which is made primarily from sodium carboxymethyl cellulose.

How it works :

When Cut-eeze gauze comes in contact with the bleeding area it quickly absorbs the water contents in the exudates to form a viscous gel which fills the wound and stops the bleeding. This process includes joining with the bloods platelets and the clot to plug the wound and close off the blood vessels.

Being a natural product the gauze converts into a gel then dissolves to saline and glucose.

What are the advantages of Cut-eeze ?

- Instant reaction with blood, stops bleeding fast
- No residue in the wound site
- No secondary bleeding
- No adverse reactions
- Ease of use and application with no pain on removal
- Safe, with no side effects
- Works with people using blood thinning medications or bleeding disorders, including haemophiliacs.

• Our Cut-eeze Haemostatic spray is a wound congealant which is ideal for small cuts and grazes.

- Spray on to a bandage or dressing, and apply to the wound
- Can also be sprayed direct to the wound area
- When the wound heals, remove residual powder and dress with plaster or bandage as usual.

