



BEESWIFT
FOCUSED ON SAFETY

Product Code

CM0563

DATASHEET

CLICK MEDICAL CUT-EEZE HAEMOSTATIC GAUZE ROLL



First aid training and regulations now include the use of haemostatic dressings and tourniquets for workers who are in high risk environments where serious trauma injuries and bleeding are a possibility.

Cut-eeze is our new range of haemostatic dressings

Description : Our haemostatic(or wound clotting) dressings comprise of a white knitted fabric which is made primarily from sodium carboxymethyl cellulose.

How it works :

When Cut-eeze gauze comes in contact with the bleeding area it quickly absorbs the water contents in the exudates to form a viscous gel which fills the wound and stops the bleeding. This process includes joining with the bloods platelets and the clot to plug the wound and close off the blood vessels.

Being a natural product the gauze converts into a gel then dissolves to saline and glucose.

What are the advantages of Cut-eeze ?

- Instant reaction with blood, stops bleeding fast
- No residue in the wound site
- No secondary bleeding
- No adverse reactions
- Ease of use and application with no pain on removal
- Safe, with no side effects
- Works with people using blood thinning medications or bleeding disorders, including haemophiliacs.

This is a 3.7m x 7.5cm Z fold gauze roll used to stop bleeding in the most severe of wounds and lacerations. Used to pack and seal larger wound areas.

Can be cut to size to use on smaller wounds

Supplied in vacuum sealed foil pouch

